

# PH NEUTRALISING EYE WASH

PART NUMBER: EWPHN500



**DESCRIPTION:**

- Recommend wash time: 2 min
- 4.9 % sterile phosphate buffer for neutralisation of acid and base that come into contact with the eye
- Washing with pH Neutral is recommended to be followed by washing with MEDIQ Eyewash
- Easily manageable bottles
- Ergonomic eye-cups for effective flow
- User-friendly design allows the bottles to be quickly opened with a single hand
- Disposable bottles; do not reuse
- CE-labelled bottles, products with 3-year expiry date

**KEY INFORMATION:**

**TAMPER-RESISTANT**

- MEDIQ buffer eyewash is made by blow-fill-seal technology which features a tamper-resistant cap to ensure quality and safety of the eye wash solution and can be easily removed in an emergency.

**INGREDIENTS**

- 4.9% Phosphate

**PRODUCT DESCRIPTION**

- MEDIQ buffer eyewash contains a 4.9% sterile phosphate buffer solution, which neutralises the eye’s fluid after accidents involving acids or alkali. The bottle is easy to take with you or can be used to refill the neutralizing eyewash station.

**SAFETY INFORMATION**

- There are no known side-effects associated with the use of sterile sodium chloride solution 0.9% w/v. Any unpleasant effects should be reported to a doctor, pharmacist or nurse.

**PRODUCT OVERVIEW:**

<b>Material:</b> Medical Grade Plastic	<b>Volume:</b> 500ml	<b>Sold As:</b> Each	<b>Colour:</b> Clear	<b>Pieces:</b> 1 (Contents Below)

## CONTENTS & QUANTITY

Medical Grade Plastic Bottle, 500ml	1
Eye-cup	1

## INGREDIENTS

4.9 % Sterile Phosphate Buffer	4.9%
--------------------------------	------

\*All contents maybe subject to change